

Anti- BL-CAM (Phospho-Tyr807) Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A100296-01	50 µl
A100296-02	100 µl

Specificity	Anti- BL-CAM (Phospho-Tyr807) (human)
Source	Rabbit Polyclonal
Application	WB ELISA
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P20273

Other Names: B-cell receptor CD22; B-lymphocyte cell adhesion molecule; BL-CAM; Leu-14; Siglec-2

Specificity and Sensitivity

BL-CAM (Phospho-Tyr807) antibody detects endogenous levels of BL-CAM only when phosphorylated at tyrosine 807.

Source and Purification

The antiserum was produced against synthesized phosphopeptide derived from human BL-CAM around the phosphorylation site of tyrosine 807 (G-D-Y^P-E-N).

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 ELISA: 1:40000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of extracts from K562 cells, treated with Na₃VO₄ (0.3nM, 40mins), using BL-CAM (Phospho-Tyr807) antibody.