

Anti-Phospho-Tyrosine Hydroxylase (Ser19) Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A100540-01	50 µl
A100540-02	100 µl

Specificity	Anti- Phospho-Tyrosine Hydroxylase (Ser19) (human mouse rat)
Source	Rabbit Polyclonal
Application	WB IHC IF/ICC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P07101

Other Names: EC 1.14.16.2, TH isoform 3, TH isoform a, TH-4, TY3H, TYH, Tyrosine 3-hydroxylase, Tyrosine 3-monoxygenase, tyrosine hydroxylase

Specificity and Sensitivity

Phospho-Tyrosine Hydroxylase (Ser19) Antibody detects endogenous levels of Tyrosine Hydroxylase only when phosphorylated at Serine 19.

Source and Purification

A synthesized peptide derived from human Tyrosine Hydroxylase around the phosphorylation site of Serine 19. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:
WB: 1:500~1:3000 IHC: 1:50~1:200 IF/ICC: 1:100-1:500


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of Tyrosine Hydroxylase phosphorylation expression in PC12 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Related Products

- PW001: Super ECL Assay kit
- E030120 : HRP, Goat Anti-Rabbit IgG(H+L)
- E030220 : AP, Goat Anti-Rabbit IgG(H+L)
- E021010: Anti-GAPDH Mouse Monoclonal Antibody
- E021020: Anti-beta Actin Mouse Monoclonal Antibody
- E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP