

Anti-Phospho-NF- kappaB p105/p50 (Ser927) Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A100560-01	50 µl
A100560-02	100 µl

Specificity	Anti- Phospho-NF- kappaB p105/p50 (Ser927) (human mouse)
Source	Rabbit Polyclonal
Application	WB IHC IF/ICC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P19838

Other Names: DNA-binding factor KBF1, EBP-1, NF-kappa-B1 p84/NF-kappa-B1 p98, NFKB1, NFkB-p50, Nuclear factor NF-kappa-B p105 subunit, nuclear factor NF-kappa-B p50 subunit

Specificity and Sensitivity

Phospho-NF- kappaB p105/p50 (Ser927) Antibody detects endogenous levels of NFkappaB p105/p50 only when phosphorylated at Serine 927.

Source and Purification

A synthesized peptide derived from human NFkappaB p105/p50 around the phosphorylation site of Serine 927. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:200 IF/ICC: 1:100-1:500


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of NF- kappaB p105/p50 phosphorylation expression in LPS treated HeLa whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP