

Anti-Fibronectin 1 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300066-01	50 µl
A300066-02	100 µl

Specificity	Anti- Fibronectin 1 (human mouse rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC
Form	Liquid, 1 mg/ml

Product

Swiss-Prot No.: P02751

Other Names: CIG, DKFZp686F10164, DKFZp686H0342, DKFZp686I1370, DKFZp686O13149, FINC, MSF, LARGE, EXTERNAL, TRANSFORMATION-SENSITIVE PROTEIN, LETS FNZ, cold-insoluble globulin, fibronectin 1 isoform 1 preproprotein, fibronectin 1 isoform 2 proprotein, fibronectin 1 isoform 3 preproprotein, fibronectin 1 isoform 4 preproprotein, fibronectin 1 isoform 5 preproprotein, fibronectin 1 isoform 6 preproprotein, fibronectin 1 isoform 7 preproprotein, migration-stimulating factor

Specificity and Sensitivity

Fibronectin 1 antibody detects endogenous levels of total Fibronectin 1 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from C-terminal of human Fibronectin 1. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:500~1:3000 IHC: 1:50~1:100 ELISA: 1:1000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions


Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human thyroid gland tissue, using Fibronectin 1 antibody.



Peptide - +

Western blot analysis of extracts from HeLa cells, using Fibronectin 1 antibody.

Related Products

- PW001: Super ECL Assay kit
- E030120 : HRP, Goat Anti-Rabbit IgG(H+L)
- E032221: Dylight 488, Donkey Anti-Rabbit IgG(H+L)
- E021010: Anti-GAPDH Mouse Monoclonal Antibody
- E021020: Anti-beta Actin Mouse Monoclonal Antibody
- E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP