

Anti- HDAC5 Polyclonal Antibody



<u>Catalog No.</u>	<u>Size</u>
A300094-01	50 µl
A300094-02	100 µl

Specificity	Anti- HDAC5 (human Mouse Rat)
Source	Rabbit Polyclonal
Application	WB ELISA IHC IF
Form	Liquid, 1 mg/ml

Specificity and Sensitivity

Swiss-Prot No.: Q9UQL6

All Names: Histone deacetylase 5, HD5, Antigen NY-CO-9

Specificity and Sensitivity

HDAC5 antibody detects endogenous levels of total HDAC5 protein.

Source and Purification

The antiserum was produced against synthesized peptide derived from human HDAC5.

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:1000~ 1:3000 IHC: 1:50~1:100

IF: 1:100~1:500 ELISA: 1:5000


Storage Buffer

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of extracts from HepG2 cells, treated with PMA (125ng/ml, 30mins), using HDAC5 antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue, using HDAC5 antibody.