

## 100 bp plus DNA Ladder



### Ordering Information

Catalog No.	100 bp plus DNA Ladder
BM002-01	50 T/250 $\mu$ l
BM002-02	4 $\times$ 50 T/250 $\mu$ l

### Description

This product is composed of 14 linear double-stranded DNA bands, 100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 1500, 2000, 3000, 5000bp. The concentration of the 500 bp band is 20 ng/ $\mu$ l, others are 10 ng/ $\mu$ l, this can be used for the quantitative analysis of DNA fragments.

### Storage Instructions

Stable for 6 months at 4  $^{\circ}$ C, -20  $^{\circ}$ C for long-term preservation.

### Images

