

D2000 DNA Marker



Ordering Information

Catalog No.	D2000 DNA Marker
BM003-01	50 T/250 μ l
BM003-02	4 \times 50 T/250 μ l

Description

This product is composed of 6 linear double-stranded DNA bands, 100, 250, 500, 750, 1000, 2000 bp. The concentration of the 750 bp band is 20 ng/ μ l, others are 10 ng/ μ l, this can be used for the quantitative analysis of DNA fragments.

Storage Instructions

Stable for 6 months at 4 $^{\circ}$ C, -20 $^{\circ}$ C for long-term preservation.

Images

