

Ordering Information

Catalog No.	DNA Marker I
BM005-01	50 T/250 μ l
BM005-02	4 \times 50 T/250 μ l

Description

This product is composed of 6 linear double-stranded DNA bands, 100, 200, 300, 400, 500, 600 bp. The concentration of the 400 bp band is 20 ng/ μ l, others are 10 ng/ μ l, this can be used for the quantitative analysis of DNA fragments.

Storage Instructions

Stable for 6 months at 4 $^{\circ}$ C, -20 $^{\circ}$ C for long-term preservation.

Images

