

## DNA Marker II



### Ordering Information

Catalog No.	DNA Marker II
BM006-01	50 T/250 $\mu$ l
BM006-02	4 $\times$ 50 T/250 $\mu$ l

### Description

This product is composed of 6 linear double-stranded DNA bands, 100, 300, 500, 700, 900, 1200 bp. The concentration of the 700 bp band is 20 ng/ $\mu$ l, others are 10 ng/ $\mu$ l, this can be used for the quantitative analysis of DNA fragments.

### Storage Instructions

Stable for 6 months at 4 °C, -20 °C for long-term preservation.

### Images

