

Ordering Information

Catalog No.	DNA Marker III
BM007-01	50 T/250 μ l
BM007-02	4 \times 50 T/250 μ l

Description

This product is composed of 7 linear double-stranded DNA bands, 200, 500, 800, 1200, 2000, 3000, 4500 bp. The concentration of the 1200 bp band is 20 ng/ μ l, others are 10 ng/ μ l, this can be used for the quantitative analysis of DNA fragments.

Storage Instructions

Stable for 6 months at 4 °C, -20 °C for long-term preservation.

Images

