

## Ordering Information

Catalog No.	DNA Marker V
BM008-01	50 T/250 $\mu$ l
BM008-02	4 $\times$ 50 T/250 $\mu$ l

## Description

This product is composed of 6 linear double-stranded DNA bands, 200, 400, 700, 1000, 1500, 2000 bp. The concentration of the 1000 bp band is 20 ng/ $\mu$ l, others are 10 ng/ $\mu$ l, this can be used for the quantitative analysis of DNA fragments.

## Storage Instructions

Stable for 6 months at 4  $^{\circ}$ C, -20  $^{\circ}$ C for long-term preservation.

## Images

