

Anti-HA Tag Mouse Monoclonal Antibody (Mix)



<u>Catalog No.</u>	<u>Size</u>
E022011-01	100µl
E022011-02	500µl
E022011-03	50µl

Specificity	Anti- HA Tag Mouse Monoclonal Antibody (Mix)
Source	Mouse Monoclonal
Application	WB IP IF
Form	Liquid, 1 mg/ml

Product

Other Names: HA epitope tag, HA1, HA2, Hemagglutinin

Specificity and Sensitivity

The HA tag antibody can recognize C-terminal, internal, and N-terminal HA-tag fusion proteins.

Source and Purification

The antiserum was produced against synthesized HA tag. The antibody was affinity-purified by affinity-chromatography using epitope-specific immunogen.

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows:

WB: 1:5000~1:10000 IP: 1:200 IF: 1:200

Storage Buffer

Mouse IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of His-MBP-HA recombinant protein with HA Tag Mouse Monoclonal Antibody (Mix) diluted at 1:10,000.

Related Products

PW001: Super ECL Assay kit

E030120 : HRP, Goat Anti-Rabbit IgG(H+L)

E030220 : AP, Goat Anti-Rabbit IgG(H+L)

E021010: Anti-GAPDH Mouse Monoclonal Antibody

E021020: Anti-beta Actin Mouse Monoclonal Antibody

E022330: Anti-His Tag Mouse Monoclonal Antibody-HRP