# **Anti-E2 Tag Monoclonal Antibody**

 Catalog No.
 Size

 E022070-01
 100μl

 E022070-02
 500μl

 E022070-03
 50μl



**Specificity** Anti- E2 Tag (C-terminal internal and N-terminal)

Source Mouse Monoclonal

Application WB

Form Liquid, 1 mg/ml

### **Application Notes**

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: Western blot (1:1000-1:10,000),

### **Source and Purification**

This monoclonal antibody is produced by immunizing mice with a synthetic peptide GVSSTSSDFRDR coupled to KLH. Antibodies are purified by protein A affinity chromatography.

### **Specificity and Sensitivity**

Anti-E2 Tag Monoclonal Antibody recognizes C-terminal, internal and N-terminal E2-tagged proteins.


#### **Storage Buffer**

PBS, pH 7.4 with 0.05% sodium azide, 50% Glycerol.

## **Storage Instructions**

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

### **Images**



Lane 1: GST-E2 fusion protein (27KD); MBP-E2 fusion protein (41KD)

Lane 2: GST protein; MBP protein

