Anti-RFP Tag Monoclonal Antibody

 Catalog No.
 Size

 E022120-01
 100μl

 E022120-02
 500μl

 E022120-03
 50μl



Specificity Anti- RFP Tag (C-terminal and N-terminal)

Source Mouse Monoclonal

Application WB

Form Liquid, 1 mg/ml

Application Notes

Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. Suggested starting dilutions are as follows: Western blot (1:1000-1:10,000),

Source and Purification

This Monoclonal antibody is produced by immunizing mice with purified RFP expressed in E.coli. Antibodies are purified by protein A affinity chromatography.

Storage Buffer

PBS, pH 7.4 with 0.05% sodium azide, 50% Glycerol.

Storage Instructions

Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.

Images



Western blot analysis of RFP-tag mouse mAb diluted at 1:5000

