

HRP Conjugated Streptavidin



<u>Catalog No.</u>	<u>Size</u>
E030100-01	100µl
E030100-02	500µl

Product Type	Secondary Antibody Related
Physical State	Liquid
Suggested Dilution Range	1:500 - 1:1000 for immunohisto/cytochemistry 1:500 - 1:1000 for ELISA and Western blotting with chromogenic substrates 1:1000 - 1:10000 for Western blotting with ECL substrates.
Storage Buffer	0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH7.6
Stabilizer	15 mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)
Preservative	None (Warning: Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.)
Storage Instructions	Store at -20°C. Avoid repeated freeze / thaw cycles.
Note	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Nothing disclosed herein is to be construed as a recommendation to use this product in violation of any patents. The information presented above is believed to be accurate. However, said information and product are offered without warranty or guarantee since the ultimate conditions of use and the variability of the materials treated are beyond our control. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.