

Revelation III 医用显微镜



目镜头部

- 双目目镜 (有三目观察头)
- 目镜倾角 30°, 可360°旋转
- 屈光度调节
- 10X/18 WF 目镜(10X/20X Super WF HP 目镜)
- 可用单目, 三日, 双目及45°滑动头
- 瞳距范围:
 - Binoc 目镜: 48-75mm
 - Trinoc 目镜: 55-75mm

物镜转换器

- 配有四个物镜适配位置
- 滚珠轴承安装 (物镜可以拧上物镜转换器上)

物镜

- DIN 消色差物镜
- 30 抗菌保证 (抗菌材质)
- 4X, 10X, 40X, 100XR (油镜)
- 自动齐焦及彩色编码

载物台

- 机械载物台平台 (140mm x 140mm)
- 同轴驱动控制
- 载物台移动范围: 73mm x 43mm
- 载物台锁定拉杆
- 缓闭液压滑动手指

聚焦

- 粗调: 30mm
- 微调: 可精确到 2µm
- 张力控制旋钮

照明

- 可以动的Abbe 聚光器, NA 1.25,
- 光圈
- 100-240V AC / 50-60Hz 自动切换输入

光源可选:

- 可变卤素光源 (12V/20W灯泡) - 可选
- 可变的LED光源 (3W球泡灯) - 标准

Revelation III

Model #	Description
R3M-BN4A-DAL3	Binoc head, 4-10-40-100 DIN Achromatic objectives, LED illumination
R3M-TN4A-DAL3	Trinoc head, 4-10-40-100 DIN Achromatic objectives, LED illumination
R3M-BN4A-DPL3	Binoc head, 4-10-40-100 DIN Plan objectives, LED illumination
R3M-TN4A-DPL3	Trinoc head, 4-10-40-100 DIN Plan objectives, LED illumination
R3M-MN4A-DAL3	Monoc head, 4-10-40-100 DIN Achromatic objectives, LED illumination

- 医疗级显微镜
- 所有镜头都是抗菌材质
- 户外便携式配件

尺寸和重量

- 高度: 15" (380mm)
- 长度: 9" (230mm)
- 宽度: 7" (178mm)
- 重量: 14 lbs. (6.4kg)



P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com