

D-LCD 放大视频采集器



特点

变焦范围: 0.7X-4.5X

缩放比例: 1: 6.5

总放大倍率: 2.5X-152.4X (8" LCD)
3.2X-198.1X (10.4" LCD)

视野: 1.3mm-81.6mm

工作距离: 46mm-331mm

分辨率: 640×480 (8" LCD)
1024×768 (10.4" LCD)

视频角度: 150°

物镜: 0.3X, 0.5X, 0.75X, 1X(标准),
1.5X, 2X

电源: AC 85-265V, 50/60Hz

LED环形灯: DC 12V, 亮度可调

3D视角附件: 35°倾角 0.75X 物镜

激光指示器: 650nm

视频输出: BNC连接器/S-视频连接器

机身尺寸: 76mm

可通过数据线连接视屏采集器和显微镜

D-LCD Zoom

Model #	Description
VSS-LC08-TS71	LCD Video Microscope (8" LCD screen)
VSS-LC10-TS71	LCD Video Microscope (10" LCD screen)



P 770.270.1394

F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com